IVD

CE

3diag - C1In - 800

REF TD-42581

LOT **581-15** 2026-11

(EN) C1 (Esterase) Inhibitor - for *IMMAGE*® 800
INSTRUCTIONS FOR USE

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<https://www.3diag.com/005>

for IFU, Scan or follow link,
and select LOT

IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA

TRIMERO Diagnostics, SL

c. València 558, 4t 2a - 08026 Barcelona (Spain)
☎ +34 93 244 86 79 - www.3diag.com

INSTRUCTIONS FOR USE

Reagents for professional use,
for *In Vitro* use only in clinical laboratory (IVD)

3diag - C1In - 800

**C1 (Esterase) Inhibitor
for IMAGE® 800**

REF TD-42581

(Product included in **REF** TD-42580)

INTENDED USE

Quantitative determination of Complement C1 (Esterase) Inhibitor (C1In, CEI), in human serum, on IMAGE® 800 Immunochemical Systems of Beckman Coulter (IMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA).

PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

CONTENTS - COMPOSITION - PREPARATION

- UDR IMAGE® Reagent Cartridge: **REAG UDR C1In**
REF TD-42581-UDR ▼ 100 test containing:
 - Antiserum Reagent: **CONT** 2.5 ml (Compartment A)
Solution of anti-human C1In antibodies.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN₃).

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. When using the reagents, it is advisable to use anti-evaporation rubber stoppers (**REF** 447169, from Beckman Coulter).
The indicated shelf life must be taken as a guideline given that, obviously, it depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

MATERIALS NEEDED, NOT SUPPLIED

- Beckman Coulter's IMAGE® 800 Immunochemical System, and accessories: dilution segments, cuvettes, etc..
- BUF1 IMAGE®, from Beckman Coulter **REF** 447650
- DIL1 IMAGE®, from Beckman Coulter **REF** 447640
- 3diag - C1In - Cal Set **REF** TD-42582
- 3diag - C1In - Control **REF** TD-42583

SAMPLES

Fresh Serum.

Samples with presence of fibrin should be centrifuged.

Do not use hemolyzed, lipemic or contaminated samples.

In bibliography⁽¹⁾ it is reported the following stability in serum:

Refrigerated: 72 hours - Frozen: 14 days

PROCEDURE

Follow the instructions⁽²⁾ for use of the analyzer to program and calibrate UDR (User Defined Reagent) assays, with the recommended parameters that are detailed below.

Assay Parameters

- Units: mg/dl
- Protocol: Non-competitive Nephelometry
- Sample or Dilution Volume: 20 µl
- Reaction Buffer Volume: 195 µl (use BUF1)
- Compartment A Volume: 20 µl
- Compartment B Volume: 0 µl (not used)
- Gain: 3
- Calibration Dilution: 1:7 (use DIL1)
- Sample Dilution: 1:7 (use DIL1)
- Reaction Time: 2.5 minutes

Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with DIL1, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range.

Calibration Parameters

- Calibrators: Use the 3diag - C1In - Cal Set
- Levels: 6
- Replicates: 2
- Model: Logistic 4 Parameters

The assays must be recalibrated, at least when a new batch of reagents is used.

The use of the calibration update option with a single calibrator level is discouraged.

PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay are given in the Technical Report, available on the website (www.3diag.com) or upon request to the Customer Support Service (✉ support@3diag.com - ☎ +34 93 244 86 79).

QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of 3diag - C1In - CONTROL.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

TRACEABILITY

Given that certified reference materials are not available, the values are referred to internal standards of highly purified proteins. Traceability is ensured by measuring the C1In in the *International Ref. Preparation for human serum complement Factors (NIBSC code: W1032)* of the *WHO (World Health Organization)*.


REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography reports, depending on the method used and the population analyzed, reference values of between: 19 to 37 mg/dl⁽¹⁾, 15 to 35 mg/dl⁽³⁾ and 17.4 to 24.0 mg/dl⁽⁴⁾.

Analyzing serum samples of 179 patients from Spanish hospital laboratories, without discarding those that presented alterations of other components of the complement, the following results have been obtained (in mg/dl):

mean	SD	range	95 percentile	90 percentile
20.7	3.91	9.18 - 38.5	14.2 - 28.5	15.0 - 26.8

In view of the results, a concentration lower than about 15 mg/dl can be taken as a significant value, indicative of a deficiency, since high C1In values do not have an established clinical significance⁽¹⁾⁽⁴⁾.

CLINICAL SIGNIFICANCE

C1 inhibitor deficiency is associated with hereditary or acquired angioedema. In 85% of cases their levels are decreased (5-30% of their normal value), while in the remaining 15%, the levels are normal but the protein is not functional. The additional measurement of C1q is key to distinguish between hereditary or acquired angioedema, since it presents normal levels in the hereditary form and low levels in the acquired form.

Low levels of C1 inhibitor also predispose to autoimmune diseases, especially Lupus Erythematosus (SLE), due to its effect of consumption of C3 and C4.

SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed⁽⁵⁾ by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

REAG

Reagent

UDR

UDR (User Defined Reagent)

C1In

C1 (Esterase) Inhibitor

CONT

Contents

BIBLIOGRAPHY

- (1) Mayo Medical Laboratories website (www.mayomedicallaboratories.com), date of consultation: 2nd August 2017.
- (2) "IMMAGE® Immunochemical System user's manual". IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA.
- (3) Putnam F.W. - "The Plasma Protein" - II Edition.
- (4) "Nomenclator de Laboratorio Clínico AEFA/AEBM" - ISBN: 84-486-0117-3.
- (5) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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IMMAGE® is a registered trademark of Beckman Coulter, Inc, Fullerton, CA.
 BN™ and Atellica® are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.
 Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.

CONTROLS INCLUDED IN  **Accuracy 365** Online Quality Control Further information and registration at www.accuracy365.com

Accuracy 365 aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.